

## DISPOSABLE BOUFFANT HOOD



**SHIELD<sup>me</sup>™ NON-STERILE DISPOSABLE BOUFFANT HOOD** The single use SHIELDme Non-Sterile Disposable Hood has design to ensure a perfect fit. Made of 100% polypropylene breathable SMS fabric, which gives the wearer cloth-like comfort. The hood features an extra-long yoke for maximum coverage when worn in conjunction with a coverall with collar, and features an elasticated face-opening with reinforced edges to avoid contamination entering the controlled environment.



### SPECIFICATIONS:

- Weight 25-50 g/m<sup>2</sup>
- Single elasticated front
- V neck and bib
- Light and breathable
- White colour
- Uni-Size
- Comfort fittings

### COMFORT FITTINGS

- Resistant to blood and body Fluid Penetration
- Resistant to viral penetration
- Resistant to Blood Borne Pathogens
- Resistant to Chemical, oil, bleach and 50% Sodium Hydroxide

*\*We reserve the right to change specifications of our product without prior notice.*